

# polyA RNA Isolation from *S. cerevisiae*

DeRisi Lab, UC San Francisco, June 2001

## Materials

2 ml size Sarstedt screw cap tubes  
BIO-RAD Disposable Chromatography Column (BIO-RAD cat. #732-6008)  
Oligo dT cellulose  
RNase-free water  
5.0 M NaCl  
0.5 M EDTA  
1.0 M Tris, pH 7.4  
10% SDS  
70% EtOH

## Before starting protocol

1. Measure out 75 mg of oligo dT cellulose into a 2 ml Sarstedt screw cap tubes.
2. Make fresh buffers:

	1X NETS	2X NETS	1X ETS
RNAse free water	10.4 ml	1.338 ml	2.85 ml
5.0 M NaCl	1.5 ml	480 $\mu$ l	-
0.5 M EDTA	250 $\mu$ l	80 $\mu$ l	60 $\mu$ l
1.0 M Tris, pH 7.4	125 $\mu$ l	40 $\mu$ l	30 $\mu$ l
10% SDS	250 $\mu$ l	62.5 $\mu$ l	60 $\mu$ l

**Note:** Amounts are for 2 preps; scale up as necessary.

- 2X NETS may need to be heated to go into solution.
- 1X ETS should be heated to 65 °C for use.

## Protocol

1. Resuspend at least 2 mg (up to 4mg) of total RNA in 750  $\mu$ l RNase-free water. In general, inputting more total RNA will yield more polyA RNA. This protocol is most successful with inputs of at least 2 mgs. (Less than 1mg of total RNA input is NOT recommended for this protocol.)
2. Wash the oligo dT cellulose 3X with 750  $\mu$ l 1X NETS. (Spin down gently between washes and aspirate off the wash.)
3. Aspirate off final 1X NETS wash and resuspend the cellulose in 750  $\mu$ l 2X NETS.
4. Incubate the RNA at 65 °C for 10 minutes.

5. Add the RNA to the cellulose and mix on rotator for 1 hour at room temperature.
6. Prepare a BIO-RAD Disposable Chromatography Column (BIO-RAD cat. #732-6008) by washing it out once with 750  $\mu$ l 1X NETS.
7. Gently pour the RNA/cellulose mixture into the column and allow it to settle by gravity. (DO NOT MIX)
8. Gently wash the column 3X with 750  $\mu$ l 1X NETS. (Pipette down the side of the column without disrupting the cellulose.)
9. Elute the mRNA into a fresh Eppendorf tube by pipetting 650  $\mu$ l of 65  $^{\circ}$ C 1X ETS directly and forcibly into the oligo dT cellulose. Repeat 1X into a second tube.
10. Add 650  $\mu$ l 3M NaAcetate to all samples. Fill the remainder of the tube with room temperature isopropanol ( $\sim$ 700  $\mu$ l). Mix well by inverting several times.
11. The sample can now sit at  $-20^{\circ}$  C to precipitate as long as desired (up to overnight). To proceed, spin sample in a 4  $^{\circ}$ C microfuge at max speed for 30 minutes to 1 hour.
12. Aspirate off supernatant and wash pellet with 250  $\mu$ l 70% EtOH and spin down 5 s.

## **polyA RNA Isolation**

13. Aspirate off the EtOH and let pellet air dry (or speed vac dry on low heat).
14. Resuspend in 15-20  $\mu$ l RNase-free water.
15. Spin down for 30 seconds and quantitate a 1:100 dilution at 260nm.

**Note:** Typically yields should be no greater than 2% of the total RNA input. Each microarray requires at least 1.5  $\mu$ g of polyA RNA/channel. Inputs less than this may lead to poor signal and additional artifacts that may contaminate datasets. For organisms other than *S. cerevisiae*, it is recommended that users optimize input RNA amounts by hybridizing the same sample at different inputs against itself: i.e., sample A 2  $\mu$ g vs. sample A 1  $\mu$ g. Optimize input amounts based on adequate signal and mean ratio=1.0.

16. Store mRNA at  $-80^{\circ}$  C.
17. Before using for reverse transcription and amino-allyl coupling, spin the sample for 30 seconds to pellet any residual oligo dT cellulose in the tube.